## DFM19HBV-SDI DFM19HBV-HDI

### USERS MANUAL

Medical Grade 19" Endoscopy Display

## Table of Contents

Table of Contents .....	2
Illustrations.....	2
About This Manual.....	3
Application and Functional Description .....	4
Some Basic Facts.....	4
Intended Use .....	4
Functional Description .....	4
Product Features DFM19HBV-SDI .....	4
Product Features DFM19HBV-HDI .....	4
Explanation of Symbols Used On The Devices, Safety Information .....	5
Explanation of information and symbols on the monitor nameplate. ....	5
General Notes on Electrical Safety .....	6
To connect the cables .....	11
To connect the cables .....	12
Cable Management (only with stand OP-DFM19-ENDO).....	13
Kensington® Anti-Theft Security Lock Slot .....	13
Adjust the Two-Way Ergonomic Stand (only with stand OP-DFM19-ENDO) .....	14
Controls & Functions .....	15
The OSD Wheel Knob .....	15
On-Screen Display .....	16
OSD – Main Menu Display .....	16
Menu Picture (VGA or DVI) .....	16
Menu Picture (Video).....	17
Menu Color Mode .....	17
Menu Source (MD-DFM19HBV-SDI) .....	17
Menu Source (MD-DFM19HBV-HDI) .....	18
Menu Management.....	18
Menu Memory Setting .....	18
Cleaning and Disinfection.....	19
Equipment Surface .....	19
Connection Cables .....	19
Appendix.....	20
Technical Information*.....	20
Specification of the medical flat screen DFM19HBV.....	20
Timings .....	21
External Dimensions.....	22
Support and Service .....	23
Repairs: .....	23
Product Warranty.....	23
Disposal.....	24

## Illustrations

Illustration 1 – Type Rate Label of the Display.....	5
Illustration 2 – Connectors and cabling of DFM19HBV-SDI.....	11
Illustration 3 – Connectors and cabling of DFM19HBV-HDI .....	12
Illustration 4 – Cable Management.....	13
Illustration 5 – Adjust the two-way ergonomic stand .....	14
Illustration 6 – OSD Wheel Knob and LED Indicator.....	15
Illustration 7 – Timings.....	21
Illustration 8 – External Dimensions (landscape front view).....	22
Illustration 9 – External Dimension (right side view) .....	22
Illustration 10 – External Dimensions (left side view).....	22
Illustration 11 – External Dimension (portrait front view).....	22

Thank you very much for purchasing the DFM19HBV, a medical safety approved TFT-LCD Display. In order to use your display correctly, please read the users manual carefully and take notice of the danger, warning and caution notices.

## About This Manual

This user manual is intended as an aid in setting up and operating the flat-screen display. We have thoroughly verified all information and hints.

The safety information given in this manual is classified as follows:

**DANGER:**

**Indicates an imminent hazard. If not avoided, the hazard will result in death or serious injury.**

**WARNING:**

**Indicates a hazard. If not avoided, the hazard can result in death or serious injury.**

**CAUTION:**

**Indicates a potential hazard. If not avoided, this hazard may result in minor personal injury or product/property damage.**

### Manufacturer

Richardson Electronics GmbH  
Division A.C.T. Kern  
Raiffeisenstr. 5  
D-78166 Donaueschingen  
Germany

Tel: +49 (0)771 / 8300-0  
Fax: +49 (0)771 / 8300-80  
Email: [medical@actkern.de](mailto:medical@actkern.de)  
Internet: [www.actkern.info](http://www.actkern.info) , [www.flatmedic.de](http://www.flatmedic.de)

Specifications are subjected to change without notice.

Revision	Document Name	Author	Date
1.2	Manual DFM19HBV-SDI V1.2_en	MK/GA	2007-10-10
1.1	Manual DFM19HBV-SDI V1.1_en	MK/GA	2007-05-30
1.0	Manual DFM19HBV-SDI V1.0_en	MK/GA	2007-01-11

## Application and Functional Description

### ***Some Basic Facts***

The medical safety approved DFM19HBV is a versatile monitor designed for ergonomic and easy handling. The liquid-crystal display (LCD) monitor supports the most common resolutions from 640x480 (VGA) up to 1280x1024 (SXGA) and presents high-brightness, sharp, low-radiation images. With the OSD wheel knob on the bottom right side of the display, the display parameters can be easily adjusted.

### ***Intended Use***

The intended use of the DFM19HBV is to display medical imaging data and endoscopy applications. It is not allowed to operate the DFM19HBV with power supplies other than the supplied LIEN Electronics AC Adapter Model LE-0312B130W or Medi Power AD-2480MB. Only this combination of display and power supply is approved for use in the medical environment. The flat-panel display and the power supply must not be used outdoors or in areas where an explosion hazard may occur. The DFM19HBV display is a medical device as defined in the Council Directive 93/42/EEC about medical devices and is class I equipment (MDD).

#### **CAUTION:**

**The VESA 100 adapter according to VESA standard is available on the rear cover for wall or swivel-arm mounting. The user is responsible for correct installation of the display.**

### ***Functional Description***

The medical flat-screen display is built into a slim and ergonomic housing. The monitor displays 16.7 Million colors with a resolution of 1280x1024 pixels. The DFM19HBV is equipped with a IPX1 protection against water drops.

The display is controlled either analog via a 15-pin D-sub VGA signal cable, digital via a DVI cable or video via different video cables. The signal is equivalent to the analog and digital standard signals of your PC. After connection of the signal the flat-screen display automatically adapts, as far as possible, to the VGA signal and presents a stable and centered image. Additional display parameters, such as brightness and contrast, can be adjusted via the on-screen menu.

#### ***Product Features DFM19HBV-SDI***

- 19" TFT Liquid Crystal Display
- protection glass, IPX1
- high brightness
- high contrast ratio
- remote control port
- power management
- 170° / 170° viewing angles
- Optimum resolution 1280 x 1024 (SXGA)
- 16.7 million colors
- VGA and DVI input
- S-Video, Composite Video, YPbPr/YCbCr, SD-SDI
- VESA 100
- user-friendly OSD-Menu with intuitive wheel-knob
- DDC2B Plug & Play
- Kensington Lock

#### ***Product Features DFM19HBV-HDI***

- 19" TFT Liquid Crystal Display
- protection glass, IPX1
- high brightness
- high contrast ratio
- remote control port
- power management
- 170° / 170° viewing angles
- Optimum resolution 1280 x 1024 (SXGA)
- 16.7 million colors
- DVI input (DVI-I, DVI-D)
- S-Video, Composite Video, YPbPr/YCbCr, HD-SDI
- VESA 100
- user-friendly OSD-Menu with intuitive wheel-knob
- DDC2B Plug & Play
- Kensington Lock

## Explanation of Symbols Used On The Devices, Safety Information

### *Explanation of information and symbols on the monitor nameplate.*

**Richardson Electronics GmbH**  
**Division A.C.T. Kern**  
[www.actkern.de](http://www.actkern.de)

**Product Name:** MD DFM19HBV-SDI

Model No.: GOD971FMII

Date of Mfg.:  
Serial No.:

Power Rating: 24Vdc, 3.0A

Caution: Only use Adapter  
Input AC100-240V, 50/60Hz, 1.7A,  
Lien Elec.Type LE-0312B130W

**CAUTION:** Do not open this cover. Trained service personnel only.  
Only use the adapter which was included in the delivery of this device or listed in the user's manual. Use of another adapter may cause a malfunction or danger.

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions:  
(1) this device may not cause harmful Interference, and  
(2) this device must accept any interference received, including interference that may cause undesired operation.

  
**CLASSIFIED**  
**UL**  
**US**  
**MEDICAL**  
**49CG**  
**E225344**

MEDICAL ELECTRICAL EQUIPMENT  
WITH RESPECT TO ELECTRIC SHOCK, FIRE AND  
MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH  
UL60601-1 AND CAN/CSA C22.2 NO.601.1

Classification: Class I  
No applied parts  
No AP/APG

  
**WEEE**  
**Directive**  
**2002/96/EC**

  
**IEC60601-1**  
**EN60601-1**  
**EN60601-1-2**

  
**Tested To Comply**  
**With FCC Standards**  
**FOR HOME OR OFFICE USE**

DI MDL4XLB160
Made In Taiwan.

**Illustration 1 – Type Rate Label of the Display**

Manufacturer:	Richardson Electronics GmbH Division A.C.T. Kern Raiffeisenstraße 5 D-78166 Donaueschingen Tel.: +49 (0) 771 / 83 00-0 Fax: +49 (0) 771 / 83 00-80
Product Name: Model No:	MD-DFM19HBV-SDI / MD-DFM19HBV-HDI GOD971FMII / GOK971FMII
DFM19HBV-SDI DFM19HBV-HDI	Standard medical display + Video (SD-SDI) Standard medical display + Video (HD-SDI)
Power Rating 24Vdc 3.0A	The medical-grade flat-screen display must only be operated with the power supply Lien LE-0312B130W or Medi Power AD-2480MB that is shipped with the device and provides the required voltage. (see also item 2)
Serial No.:	Serial Number
EN 60601-1 UL 60601-1 CAN/CSA 22.2	The device fulfills the requirements of the medical standard EN 60601-1 and cUL 60601-1 and CAN/CSA 22.2 No.601.1
	The symbol indicates compliance with the medical standard EN 60601-1 and cUL 60601
	Consult User's Manual! Refer to the user manual and observe the safety information and instructions explaining how to put the monitor into service.

### **General Notes on Electrical Safety**

Connecting medical electrical equipment to other medical electrical devices and/or other non medical electrical devices, be it for diagnostic or therapeutic purposes or for monitoring, creates a medical electrical system. The standard EN 60601-1-1 and cUL 60601 describes the requirements for these systems.

These systems must be composed of devices that comply either with the standard EN 60601-1 and or cUL 60601 or with the applicable IEC / ISO standard (PCs, for instance, must meet the requirements of the standard IEC 60950) or that are equipped with additional safety features, such as a protective earth conductor, isolating transformer or electrically isolated interface (e.g. optocoupler). This applies to devices of protection class I or II. It is not permitted to integrate protection class 0 devices in medical electrical systems.

In the patient environment, the following maximum enclosure leakage currents (current flowing between two device enclosures or between the enclosure of a device and the protective earth conductor) are permitted for medical devices:

normal condition: 0.1 mA

single fault condition: 0.5 mA (e.g. protective earth conductor interrupted)

Before putting the system into service for the first time, a specialist is required to check that these values are not exceeded.

If the enclosure leakage current of a device or of the system exceeds the limit values stated above, additional protective measures, such as connection to protective earth conductor, isolating transformers or conductive castors, must be implemented to correct this problem.

### **Safety Information**

#### **DANGER:**

**Explosion Hazard – The display including the power supply is not designed for use in areas of medical locations where an explosion hazard may occur.**

**In appropriate usage of flammable anesthetics, skin cleansing agents and disinfectants the medical Displays of the DFM series can be used risk less in surgical rooms and intensive care units.**

**It is not permitted to operate the devices in the presence of flammable substances (e.g. gas).**

#### **WARNING:**

**Shock Hazard — Use only following medical power supply:**

**Medical Power Supply Model Lien LE-0312B130W AC-Input: 100 – 240V~, 50 – 60Hz ,**

**DC-Output: +24V= / 5.41A**

**or Medical Power Supply Model Medi Power AD-2480MB, AC-Input: 100 – 240V~, 50 – 60Hz ,**

**DC-Output: +24V= / 3.33A**

**WARNING:**

**Personal Injury —** Medical electrical devices may be connected to other devices or to parts of systems only when it has been made certain that there is no danger to the patient, the user, or the environment as a result.

In those instances where there is any element of doubt concerning the safety of connected devices, the user must contact the manufacturers concerned or other informed experts as to whether there is any possible danger to the patient, the operator, or the environment as a result of the proposed combination of devices.

Connecting medical electrical devices to other medical electrical devices and/or other non-medical electrical devices creates a medical system which must comply with the system standard

EN 60601-1-1 / cUL 60601. Before putting the system into service, a specialist is required to ascertain its conformity with the system standard.

**WARNING:**

**Personal Injury —** Before putting the device into service, check all cables for signs of damage. Damaged cables must be replaced immediately with the original manufacturer's cables. The safety information and certificates do not apply when cables from other manufacturers are used.

Do not open the flat-screen display or perform any service interventions. This may render your warranty claims null and void. Only authorized technicians are permitted to service the flat-screen display.

Do not put the flat-screen display into service

- when liquids have entered the device,
- when it was dropped and the housing is damaged.

The flat-screen display including the power supply has no special protection against the ingress of water (dripping, splashing water). Do not immerse the devices into liquids

- Disconnect the power supply from the power line before cleaning. Use a moist cloth.
- Disconnect the flat-screen display from the power supply before cleaning.

**WARNING:**

**Equipment Failure —** Magnetic and electrical fields are capable of interfering with the proper performance of the devices. For this reason make sure that all devices operated in the vicinity comply with the relevant EMC requirements. X-ray equipment, MRI devices, radio systems, and cellular telephones are possible sources of interference as they may emit higher levels of electromagnetic radiation.

Keep the system away from these devices and verify its performance before use.

**WARNING:**

**Risk of Suffocation —** Dispose of the packaging material, observing the applicable waste-control regulations and keeping it out of children's reach.

**CAUTION:**

**Personal Injury —** All devices used must comply with the applicable IEC/EN or ISO standard (PCs, for instance, must comply with IEC 60950). This applies to custom designs in particular

**CAUTION:**

**Personal Injury, Equipment Damage —** Check that the ambient conditions specified by the manufacturers of the different devices are ensured.

Before cleaning the flat-screen display, turn it off and remove the 24VDC connector or the power supply connection cable. Do not spray liquid cleaning agents directly onto the display. Spray the cleaning agent onto a cloth and wipe the display clean without exerting any pressure.

Do not place any objects on the signal cables. Route the signal cables and the power cable to the flat-screen display such that they do not present a stumbling hazard.

**CAUTION:**

**Equipment Damage —** Do not set up the flat-screen display in the immediate vicinity of a window. Rain, water, humidity and sunlight may damage the flat-screen display.

Choose a flat, solid surface to set up the flat-screen display.

Do not exert pressure on the LCD. Excessive pressure may permanently damage the display.

For safe storage of the flat-screen display, a temperature between -20°C and +60°C must be ensured. Temperatures above or below this range may permanently damage the flat-screen display. Always store your flat-screen display in its original shipping box.

**WARNING:**

**Personal Injury —** During intracardiac application, it is recommended to connect the individual devices to the room's potential equalization system via separate potential equalization cables.

For intracardiac application, the application guidelines for medical electrical devices set forth in the German standard VDE 0753, part 2, must be observed.

## **Manufacturer's declaration of compliance with standards and directives**

### **Manufacturer**

Richardson Electronics GmbH

Division A.C.T. Kern

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D-78166 Donaueschingen

Tel.: +49 (0) 771 / 8300-0

Fax: +49 (0) 771 / 8300-80

### Manufacturer declaration of compliance with standards and regulations

Certification that the medical Monitor

#### **DFM19HBV-SDI**

#### **DFM19HBV-HDI**

is developed, manufactured and tested according to the following guidelines, which are fixed by the council of European Community:

EN 60601-1-2: 1993    EMC – Emission & EMC – Immunity

EMC-Tests have been carried out using standard cables with cables length of up to 2,5 meters. If longer cables are used, additional tests might be necessary to prove compliance with the EMC-directive.

IEC 601-1

IEC 601-1: 1988 + A1: 1991 + A2: 1995 = EN 60601-1: 1990 + A1: 1993 + A2: 1995

MEDICAL ELECTRONIC EQUIPMENT

PART1: GENERAL REQUIREMENTS FOR SAFETY

**93/42/EWG**

## Your Medical Display DFM19HBV

Before you start to install your medical display, make sure the carton box is not damaged and that all equipment is delivered completely.

### ***Transport Damage***

If the equipment has obviously been damaged in transit, make a written report immediately upon arrival of the goods and notify the shipping agent. The shipping agent must be informed of hidden damage in transit in writing within 2 work days.

### ***Unpacking***

Before unpacking the flat-screen display, prepare a suitable area for its set-up. This should be a flat, solid and clean surface in the vicinity of a wall outlet. Although the flat-screen display consumes only little energy, you should choose a location that affords sufficient ventilation. Furthermore, the air vents should not be covered to prevent overheating. Choose a location where the display will not be frequently exposed to direct sunlight or other strong sources of light. The light reflects on the display, dazzling you, and you will have problems seeing the displayed information.

After unpacking the display, check whether the shipment is complete and all items are intact: flat-screen display DFM19HBV

- User's Manual
- medical-grade power supply
- VGA (D-Sub 15 pin) standard cable
- DVI-D (DVI-D 18 pin) standard cable
- S-Video cable
- Power cord


If items are missing or damaged, contact your dealer immediately. Keep the packaging material and the shipping carton for storage or shipment of the display.

### ***Setup and positioning***

No tools are required to set up the DFM19HBV. Simply follow the instructions outlined in the next few pages. Connectors for the signal cables and power are located on the back of the panel. Please refer to the diagram on this page for the connector configuration.

Set up the DFM19HBV in a location where all connections are easily accessible.

The base of the flat-screen display requires a flat and solid surface (horizontal max. +5°) with good ventilation. Choose a location where sunlight or other bright light does not directly fall on the display. The flat-screen display should not be exposed to excessive heat, humidity or direct sunlight, because this could cause damage to the device.

**Illustration 2 – Connectors and cabling of DFM19HBV-SDI**

### To connect the cables

#### 1 24V DC connector

The power input on your flat screen is marked with DC IN 24V please put the 24V connector (4-pin) into this slot. You will feel a little click if the connection is done correctly. Connect your power supply with the power outlet.

#### 2 Analog input connector (VGA)

Connect the 15pin D-Sub VGA cable to this connector and to your graphic adapter.

#### 3 Digital input connector (DVI)

Connect the DVI cable to this connector and to your graphic adapter.

#### 4 Control input connector (GPIO DIN-9)

#### 5 Component/RGBS input connector

Connect a component video cable to this connector and a component video output device (f.e. DVD, High Definition (HD), RGBS video image)

#### 6 Component/RGBS output connector

#### 7 S-Video input connector

S-Video signal input. Connect a S-Video signal cable to this connector and a S-Video output video device.

#### 8 S-Video output connector

#### 9 CVBS Video input connector

Composite signal input. Connect a composite video cable to this connector and a video output device.

#### 10 CVBS/Composite Video output connector


#### 11 SD-SDI input connector

Connect SDI video connector and a SDI video output device with a 75-Ohm coaxial cable (e.g. Camcorder)

#### 12 SD-SDI output connector

### CAUTION:

**When you disconnect the cord/cables, be sure to hold the connector and not the cable itself.**

**Illustration 3 – Connectors and cabling of DFM19HBV-HDI**

### To connect the cables

**1 24V DC connector**

The power input on your flat screen is marked with DC IN 24V please put the 24V connector (4-pin) into this slot. You will feel a little click if the connection is done correctly. Connect your power supply with the power outlet.

**2 Digital input connector (DVI-I)**

Connect the DVI cable to this connector and to your graphic adapter.

**3 Digital input connector (DVI-D)**

Connect the DVI cable to this connector and to your graphic adapter.

**4 Control input connector (GPIO DIN-9)**

**5 Component/RGBS input connector**

Connect a component video cable to this connector and a component video output device (f.e. DVD, High Definition (HD), RGBS video image)

**6 Component/RGBS output connector**

**7 S-Video input connector**

S-Video signal input. Connect a S-Video signal cable to this connector and a S-Video output video device.

**8 S-Video output connector**

**9 CVBS Video input connector**

Composite signal input. Connect a composite video cable to this connector and a video output device.

**10 CVBS/Composite Video output connector**

**11 HD-SDI input connector**

Connect SDI video connector and a SDI video output device with a 75-Ohm coaxial cable (e.g. Camcorder)


**12 HD-SDI output connector**

**CAUTION:**

**When you disconnect the cord/cables, be sure to hold the connector and not the cable itself.**

### Cable Management (only with stand OP-DFM19-ENDO)

The DFM19HBV has built-in cable guides and cover to help you organize and route the cables neatly on your desktop space. The cable cover is located on the back of the monitor stand and it is opened by gently pulling it off the stand. Under the cover, the cables can be held in position by numerous clips. There are several ways to route the cables and you may use any way you find suitable for your situation. However, the routing method illustrated by the drawing on the right is recommended to give the best result in most situations both aesthetically and functionally. Once the cables are in position, gently snap the cover back on to the stand.

#### Caution

- Do not give signal cables or power cord a strong pull.
- A strong pull may cause damage on the connection part inside the display where high-precision processing technology is applied. It also may result in felling the display and damaging on the LCD panel. In the worst-case scenario, the LCD panel will be broken.

**Illustration 4 – Cable Management**

### Kensington® Anti-Theft Security Lock Slot

The DFM19HBV is equipped with a security lock slot compatible to Kensington® security lock type. The security cable lock may be available thru your dealer or it can be purchased at most computer peripheral stores near you.

### Adjust the Two-Way Ergonomic Stand (only with stand OP-DFM19-ENDO)


Adjust the two-way (tilt and pivot) ergonomic stand of your DFM19HBV for your maximum viewing comfort. To minimize eye fatigue, we suggest you allow a minimum viewing distance of 18 inches.

**CAUTION:**

In order to protect the LCD, be sure to hold the edge of bezel whenever you adjust it and do not touch the screen.

**CAUTION:**

When setting the LCD to portrait orientation, lift the LCD panel to the top of the tilt stand as shown in the figure below and rotate it by 90 degrees.

**Illustration 5 – Adjust the two-way ergonomic stand**


- (1) Tilt the LCD panel upward
- (2) Lift the Panel to the top
- (3) Rotate the LCD panel by 90 degrees.
- (4) Replace the LCD panel to the bottom while in the portrait orientation.

## Controls & Functions

### The OSD Wheel Knob

The OSD Wheel Knob is a multi-functional device located behind the LED Indicator on the right side of the front bezel. It has three movements - rotate upward, rotate downward and press inward as a button.

- Rotate Downward : Move Up/Right, Increase
- Rotate Upward : Move Downward/Left, Decrease
- Button Press : Execute, Do, Save

**Illustration 6 – OSD Wheel Knob and LED Indicator**

- (1) Turn On and OFF  
Press the OSD wheel knob in order to turn ON the DFM19HBV if the LED indicator is off. In order to turn OFF, press the OSD wheel knob for 2 seconds until the LED indicator turns off as well.
- (2) LED Indicator  
This LED indicator turns green when the power is switched ON and the power cord is properly attached. It turns orange when the monitor goes into a power saving mode (standby). Please refer to the power management section of this manual for more information.

#### OSD - On Screen Display usage

While the monitor is turned ON (LED indicator is green and image on the screen), pressing on the Control Dial (OSD wheel knob) activates the OSD. While the OSD menu is active, use the three way movements of the control dial to adjust the settings.

- rotate downward: move up / right, increase
- rotate upward: move downward / left, decrease
- button press: Execute, Do, Save

Mode	Power consumption (only flat screen)	LED indicator
ON	Approx. 60,0 W	Green
STANDBY	Approx 15,0 W	Orange
OFF		Off, no light

## On-Screen Display

The OSD menu is activated by pressing the Control Dial inward and you can select and adjust the function of your choice by rotating and clicking the Control Dial. The main menu displays a list of submenu icons and the current video input mode. Rotate the dial to move the highlights to the control you would like to adjust, then press the Control Dial inward to select that control or to activate that function. Depending on the control you selected, a submenu of the control with a status bar will appear. Rotate the Control Dial to adjust the control.

When you have finished making the adjustments, the setting is saved automatically by exiting the control function. If you do not touch the control dial for 3 seconds, the OSD is automatically exited saving your current settings.

## OSD – Main Menu Display

The OSD menu is displayed with icons, adjustment items, and setting values as shown below. The charts below display the function tree and brief explanations of the functions.

The Main Menu consists of followings menu points:

<b>Exit</b>	Close the OSD screen
<b>Auto Setup</b>	Automatically adjusts the display according horizontal and vertical adjustments, position, frequency and phase. It takes only a couple of seconds. After that you have a perfect image. For the digital input DVI, this function is not available.
<b>Brightness</b>	Adjusts the brightness of the full screen at the range from 0 to 100. Rotating the control dial up/down to adjust the brightness any time while the OSD is off
<b>Picture</b>	The picture menu allows to make changes of the picture
<b>Color Mode</b>	Allows to select color mode or color temperature
<b>Source</b>	Allows to select the input signals of the display: VGA, DVI, S-Video, Composite, Component, RGBS and SDI
<b>Management</b>	Contains several adjustments in order to change miscellaneous properties of your display, e.g. language, OSD position etc.
<b>Memory Setting</b>	Allows to store five different user-defined settings like brightness, contrast etc.

## Menu Picture (VGA or DVI)

The picture menu allows making changes of the picture. Following sub menus are available:

<b>Exit</b>	Returns to Main Menu
<b>Contrast</b>	Adjusts the contrast of the full screen at the range from 0 to 100.
<b>Black Level</b>	Adjust black level of the screen. The larger the value is, the brighter the brightness of the background contrast is, and vice versa. Adjustable at the range from 0 to 100.
<b>Display</b>	Adjust the horizontal position of the display. When rotating Control dial clockwise, the display moves to the left, and vice versa. For the digital input DVI, this function is not available. Adjust the vertical position of the display. When rotating Control dial clockwise, the display moves upward, and vice versa. For the digital input DVI, this function is not available.
<b>Clock/Phase</b>	Make "Phase" and "Clock" adjustments, if any flickers, blurs, or horizontal stripes on screen. Adjustable at the range from 0 to 100. For the digital input DVI, this function is not available.

## Menu Picture (Video)

The picture menu allows to make changes of the picture. Following sub menus are available:

<b>Exit</b>	Returns to Main Menu
<b>Contrast</b>	Adjusts the contrast of the full screen at the range from 0 to 100.
<b>Black Level</b>	Adjust black level of the screen. The larger the value is, the brighter the brightness of the background contrast is, and vice versa. Adjustable at the range from 0 to 100.
<b>Saturation</b>	Adjust saturation of the screen. The larger the value is, the more intensive is the specific hue. Adjustable at the range from 0 to 100.
<b>Hue</b>	Adjust gradation of colour. Adjustable at the range from 0 to 100.
<b>Sharpness</b>	Adjust sharpness of the screen. The larger the value is, the outline gets sharper. To make characters look sharper, set a large value. To make pictures and images look finer, set a small value. Adjustable at the range from 0 to 10.

## Menu Color Mode

Selects the color mode or color temperature.

<b>Exit</b>	Returns to Main Menu
<b>Normal</b>	Adjust preset color temperature of Cool, Neutral, Warm or User for customized red, green and blue levels: The user can adjust and set tones. - R Adjust red and equivalent colors at the range from 0 to 100. The greater the value is, the deeper the color is, and vice versa. - G Adjust green and equivalent colors at the range from 0 to 100. The greater the value is, the deeper the color is, and vice versa. - B Adjust blue and equivalent colors at the range from 0 to 100. The greater the value is, the deeper the color is, and vice versa.
<b>Gamma</b>	Adjust of gamma value (Gamma 1.8, Gamma 2.0, Gamma 2.2, Gamma 2.4)

## Menu Source (MD-DFM19HBV-SDI)

Selects the Input signal of your display. The signal is not auto detected, so you have to choose the input signal manually between VGA, DVI, S-Video, Composite Video, Component (YPbPr/YCbCr), RGBS or SDI.

<b>Exit</b>	Returns to Main Menu
<b>Analog VGA</b>	Select D-Sub 15-pin input signal (analog)
<b>Digital</b>	Select DVI 24-pin input signal (digital)
<b>S-Video</b>	Select S-Video input signal (Video)
<b>Composite</b>	Select Composite (CVBS) input signal (Video)
<b>YPbPr/YCbCr</b>	Select Component input signal (Video)
<b>RGBS</b>	Select RGBS input signal (Video)
<b>SDI</b>	Select SDI input signal (Video)

## Menu Source (MD-DFM19HBV-HDI)

Selects the Input signal of your display. The signal is not auto detected, so you have to choose the input signal manually between DVI-I, DVI-D, S-Video, Composite Video, Component (YPbPr/YCbCr), RGBS or HD-SDI.

<b>Exit</b>	Returns to Main Menu
<b>Analog DVI-I</b>	Select DVI 24-pin input signal (analog)
<b>Digital DVI-I</b>	Select DVI 24-pin input signal (digital)
<b>Digital DVI-D</b>	Select DVI 24-pin input signal (digital)
<b>S-Video</b>	Select S-Video input signal (Video)
<b>Composite</b>	Select Composite (CVBS) input signal (Video)
<b>YPbPr/YCbCr</b>	Select Component input signal (Video)
<b>RGBS</b>	Select RGBS input signal (Video)
<b>HD/HD-SDI</b>	Select HD-SDI input signal (Video)

## Menu Management

Utility contains several adjustments in order to change miscellaneous properties of your display, e.g. language, OSD position etc.

<b>Exit</b>	Returns to Main Menu
<b>Scaling</b>	Select one scale from the followings to display images: Full, 1:1, 4:3, 16:9, Full 16:9. When using scaling function in expanding direction other than 1:1, image will be expanded and thus blur. In Full 16:9 scale, image will be expanded in height and cut in width.
<b>OSD Display</b>	Adjust the OSD position on the screen. OSD H. Position: When rotating Control dial clockwise, OSD moves to the right, and vice versa. OSD V. Position: When rotating Control dial clockwise, OSD moves downward, and vice versa.
<b>Language</b>	Selecting this control, then rotate the Control Dial to select the language you want: English, French, German, Italian, Spanish, Japanese. Press the Control Dial to execute when selected
<b>Recall</b>	Initialize the data such as display position and automatic adjustment data to factory default

## Menu Memory Setting

Allows to store five different user-defined settings like brightness, contrast etc

<b>Exit</b>	Returns to Main Menu
<b>Memory Save</b>	Saves the current settings into the display memory. You can save five different settings (USER1, USER2...)
<b>Memory Load</b>	Loads the previous user save data. Choose "default" to recall the original preset data.

## Cleaning and Disinfection

**WARNING:**

**Shock Hazard — Disconnect the power supply from the power line and all display connection cables before cleaning or disinfection.**

**Do not reconnect the devices to the power line or to other devices before all cleaned parts are completely dry.**

### Equipment Surface

- Wipe the equipment surface down with a moist cloth; do not allow liquid to enter the equipment. All cleaning agents and disinfectants commonly used in hospitals are suitable. Do not use abrasive products.

**CAUTION:**

**Personal Injury, Equipment Damage — Strictly observe the instructions for use of the cleaning agents and disinfectants supplied by the manufacturers.**

### Connection Cables

- Disconnect the cables from the device and remove the power connector before cleaning or disinfection. When disconnecting the cable, be sure to pull on the connector, not on the cable.
- **Clean** the cable by rubbing it down with a cloth moistened with soap water. Use a **disinfectant** for disinfection. Do not immerse the cable in liquid.

## Appendix

### Technical Information\*

(\* Applicable after a minimum warm-up time of 30 minutes)

#### Specification of the medical flat screen DFM19HBV

Module	19 Inch SXGA-TFT-LCD, diagonal 48,3 cm
Active Area	376 mm horizontal x 301 mm vertical
Resolution	1280 x 1024
Pixel pitch	0,294 mm x 0,294 mm
Response Time	12ms
Colors	16.7 Million
Viewing angle	178° / 178° horizontal / vertical
Contrast ratio	600:1 typical
Brightness	300cd/m <sup>2</sup> typical
Analog Input Signal	0.7 Vpp/75 Ohm Horizontal sync and composite sync signal: TTL level 2.5~5.5V (plus . minus) Vertical sync signal: TTL level 2.5~5.5V (plus . minus)
Input Signal DFM19HBV-SDI	VGA, DVI S-Video BNC connector: Composite signal BNC connector: Component signal / RGBS signal SD-SDI connector: Compliant to SMPTE-259M, NTSC, PAL
Input Signal DFM19HBV-HDI	DVI-I, DVI-D S-Video BNC connector: Composite signal BNC connector: Component signal / RGBS signal HD-SDI connector: Compliant to SMPTE-292M, NTSC, PAL
Power consumption	24V= (60W / 15W Standby)
Power management	according to VESA® DPMS
Medical Power Supply	MEDICAL POWER SUPPLY Medi Power AD-2480MB, AC-Input: 100 – 240V~, 50 – 60Hz , DC-Output: +24V= / 3.33A or Lien LE-0312B130W AC-Input: 100 – 240V~, 50 – 60Hz , DC-Output: +24V= / 5.41A
Dimension	Width 432mm x Depth 251mm x Height 466mm (landscape)
Weight	9,2 kg
Temperature Range	Operation: 5°C to ~ 35°C Storage: -20°C ~ +60°C
Humidity:	Operation: 20% ~ 80%, (no condensation) Storage: 10% ~ 90%, (no condensation)
Plug & Play	VESA® DDC2B
Approvals	UL60601-1, CSA C22.2 No. 601.1, TÜV (EN60601-1), FCC, CE
Options	Stand base OP-DFM19-ENDO

Despite high quality controls of the panel manufacturer sometimes defect pixels can not be avoided.  
Compare with specification of panel manufacturer.

## Timings

※ The display may not work correctly with timings other than listed below.

◎ : Recommended timing    ○ : Applied timing

Display mode		Resolution		Frequency		
		Horizontal	Vertical	Horizontal (kHz)	Vertical (Hz)	
VGA		720	400	31.5	70	○
		640	480	31.5	60	○
VESA		640	480	37.9	72	○
		640	480	37.5	75	○
		640	480	43.3	85	○
		800	600	35.1	56	○
		800	600	37.9	60	○
		800	600	48.1	72	○
		800	600	46.9	75	○
		800	600	53.7	85	○
		1024	768	48.4	60	○
		1024	768	56.5	70	○
		1024	768	60.0	75	○
		1024	768	68.7	85	○
		1280	1024	64.0	60	◎
		1280	1024	80.0	75	○
		1600	1200	75.0	60	△
MAC		640	480	35.0	67	○
		832	624	49.7	75	○
		1024	768	60.2	75	○
Video Input	NTSC			15.73	60	Interlace display
	PAL			15.63	50	Interlace display

### Illustration 7 – Timings

Caution: 1600x1200 is a downscaled pseudo-timing!

## External Dimensions


Illustration 8 – External Dimensions (landscape front view)


Illustration 9 – External Dimension (right side view)


Illustration 10 – External Dimensions (left side view)


Illustration 11 – External Dimension (portrait front view)

## Support and Service

Richardson Electronics GmbH, Division A.C.T. Kern from Donaueschingen/Germany is active in national and international markets. Full application support and service is available. For technical support or service problems, please contact our service team at:

Richardson Electronics GmbH  
Division A.C.T. Kern  
Raiffeisenstr. 5  
D-78166 Donaueschingen  
Germany

Tel: +49 (0)771 / 8300-0  
Fax: +49 (0)771 / 8300-80  
Email: [medical@actkern.de](mailto:medical@actkern.de)  
Internet: [www.actkern.info](http://www.actkern.info), [www.flatmedic.de](http://www.flatmedic.de)

### **Repairs:**

Equipment for repair must be returned in the original shipping carton (postage prepaid) to the above service address. RMA number and error description must be included. After receipt of the equipment for repair we will send you an order confirmation and a cost estimate (for repairs after the warranty period).

For questions about servicing and repairs, please call your contact in our Sales team. There, you will receive more detailed information about handling service and repair cases.

## Product Warranty

In case of defects, please advise Richardson Electronics GmbH according to the aforementioned RMA Handling Procedures promptly, but no later than 8 days after receipt of goods, or in case of hidden defects no later than 3 days after discovery. The customer has no right to return the goods without prior approval from Richardson Electronics GmbH.

If there has been no prior agreement, the warranty-period is 12 months commencing from the invoice date by Richardson Electronics GmbH.

Richardson Electronics GmbH will not bear costs resulting from re-performance, particularly infrastructure and transport costs, if the goods have been placed to another location than the place of delivery. The customer shall choose the lowest-price form of transportation.

Claims for warranty defects shall not exist in cases of:

- improper usage
- modification of the goods by the customer
- using the goods outside its field of application or its electrical specifications
- neglect or misuse of the goods by the customer
- natural wear and tear and damages, defects, reduced output, and changes of conditions or operation of our goods due to extraneous causes ( e.g. impact, blows, agitation, water, fire), improper storage, treatment or installation, unusual climatic conditions, special conditions at receipt or operational conditions at the location of use, or force majeure
- defects due to construction and material deficiencies, when the customer has specified the construction or the material
- image retention on the LCD-Panel, which is caused by permanent pictures as well as Gap Mura on the LCD-Panel

If there are any other issues which have not been stated, the current version of the General Terms and Conditions of Delivery and Payment issued by Richardson Electronics GmbH should apply.

## **Disposal**

For the disposal of your flat-screen display, please observe the applicable local regulations.

Richardson Electronics GmbH  
Division A.C.T. Kern  
Raiffeisenstr. 5  
D-78166 Donaueschingen  
Germany

Tel: +49 (0)771 / 8300-0  
Fax: +49 (0)771 / 8300-80  
Email: [medical@actkern.de](mailto:medical@actkern.de)  
Internet: [www.actkern.com](http://www.actkern.com)  
[www.flatmedic.com](http://www.flatmedic.com)

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© 2007 Richardson Electronics GmbH
El cliente no tendrá derecho a devolver los bienes sin una aprobación previa de Richardson Electronics GmbH.